

WHY JMEDGUARD STANDARD™?

GREATER THAN
99.99%*

kill effectiveness on E.coli, P.Aeruginosa and S.Aureus at 90 days after coating

jMedGuard™ Antimicrobial Coatings are non-toxic and durable for both home and business use.



**CONCENTRATED
SPRAY & WIPE
ANY SURFACE**



**ONE SPRAY
90 DAYS STAY**



**ALWAYS-ON
ANTIMICROBIAL**

- INVISIBLE COATING
- WATER-BASED FORMULA
- NON-CORROSIVE
- READY-TO-USE FORMULA

*Based on SGS test report following ISO 22196:2011 (measurement of antibacterial activity on plastics and other non-porous surfaces)

ABOUT JMEDGUARD STANDARD™

jMedGuard Standard™ offers the state-of-the-art technology in persistent antimicrobial coating. Distilled from years of experience in medical implant device coating from jMedTech Pte Ltd, jMedGuard Standard™ offers the new standard in durable surface self-disinfection by serving commercial and residential consumers.

OUR VISION

Guarding you. Protecting the Future.

OUR MISSION

To revolutionize the cleaning and disinfecting product sector through innovation and development of medically effective solutions for a cleaner and safer future.



jMedGuard
ANTIMICROBIAL SOLUTIONS

JMEDGUARD PTE LTD

📍 JTC Medtech Hub 08-09
2 Tukang Innovation Grove
Singapore 618305

☎ +65 6465 5909

✉ info@jmedguard.com

🌐 www.jMedGuard.com



jMedGuard
ANTIMICROBIAL SOLUTIONS

STANDARD

**SELF DISINFECTING
COATING SPRAY**



90 days protection

99.99% Kills
Bacteria & Viruses

VALUE FOR MONEY



JMEDGUARD™ STANDARD SELF-DISINFECTING COATING SPRAY (100mL)

For Travel, Home, & Office

1 Bottle Usage = **10sqm**

Area Size of **700 X**



**SAVE
MORE**



JMEDGUARD™ STANDARD SELF-DISINFECTING COATING SPRAY (300mL)

For Home & Office

1 Bottle Usage = **30sqm**

Area Size of **1000 X**



**SAVE
MORE**

JMEDGUARD™ STANDARD SELF-DISINFECTING COATING SOLUTION (5L)

For Commercial & Industrial Purposes



1 Bottle Usage = **500sqm**

Area Size of

**4 X Average
HDB FLATS**



**SAVE
ALOT
MORE**

350, 000 X



HOW TO USE JMEDGUARD™ STANDARD

1



CLEAN OILY OR SOILED SURFACES

If surface is soiled or oily,
prepare surface by
wiping it clean before
coating treatment.

2



SPRAY & WIPE

Spray a thin layer of
coating solution on the
cleaned surface.
Afterwards, wipe and
spread the coating
solution evenly with a
non-absorbent cloth and
leave the surface to dry.

3



PEACE OF MIND FOR 90 DAYS!

jMedGuard Standard
forms a long-lasting
self-disinfecting coating
upon drying and protects
against viruses, bacterias,
molds, fungi, and algae
for more than 90 days.

Regular cleaning is still necessary as coating
enhances protection in between cleaning runs,
by guarding 24/7 against recontamination as
well as cross-contaminations.


HOW JMEDGUARD™ STANDARD WORKS



High microbial count on surface



jMedGuard™ Standard is applied and
polymerises into a durable coating on surface



jMedGuard™ Standard's long molecular chains
pierces and kills germs

WHERE TO USE JMEDGUARD™ STANDARD*



jMedGuard™'s antimicrobial surface solutions
can be used in a variety of environments, from
residential use, food and beverage restaurants,
healthcare, public and private transportation,
schools, athletic equipment and other areas
with cases of high touch/contact.

*Antimicrobial surface coating durability is dependent on
method of application and physical wear on surface